

ARX Proximal Humerus System

Greater Tuberosity Humeral Fracture Plate

Indication




Designed to fix for splitting and avulsion fractures of greater tuberosity



Features

- Easy to restore the normal GT anatomy with stable fixation
- Provide to lessen deltoid damage using deltoid splitting approach
- Permits early range-of-motion exercises and preserves the integrity of the rotator cuff.
- More suture holes available
- Anatomically contoured tuberosity plating


Shape and Positioning



Types of Plates



Ordering Information



Code	Length	Width	Holes	Color	ect.
35-PLHU-A01-L	42.3 mm	24.7 mm	8	Green	Left
35-PLHU-A01-R	42.3 mm	24.7 mm	8	Blue	Right
35-PLHU-S01-L	42.3 mm	24.7 mm	8	Green	Left
35-PLHU-S01-R	42.3 mm	24.7 mm	8	Blue	Right

