


# ARX KNEE SYSTEM

## Patella Plate

### Indication

Designed to fix the fractures of the patella bone



### Features

- Anatomically pre-contoured
- More stable fixation than conventional tension band wiring
- Ø2.5/2.8 Locking and non-locking screws used
- Rigid fixation with crossed screws
- 1.3mm Thickness to minimize the foreign body sensation


Ø1.6 Guide Pin Holes

Compression Holes For Additional Reduction


Ø1.6 Guide Pin Holes

Rigid Fixation by 2.8 Cortical Screw


Crossed Screws for Rigid Fixation




## Features



## Ordering Information



Anterior Patella Plate				
Code	Type	Hole	Thickness	Length
25-ANPA-109	Small	9	1.3 mm	44.4 mm
25-ANPA-209	Large	9	1.3 mm	48.4 mm



## JEILMEDICAL