




Distal Ulna System

Rigid & Safe Anatomical Plates


Upper Extremity





Distal Ulna System

Rigid & Safe Anatomical Plates



Compatible with
Ø2.0 Locking & Non-Locking Screw



Hexalobular Internal (T7)

**Specially designed solution
for distal ulna fracture management**
Low Profile with Improved Hook Strength

Narrow and low profile design and smooth edges to minimize soft tissue irritation



Oblong hole for plate position adjustment

Compression hole for inter-fragmentary compression




1.8mm

Pointed hooks for ulnar styloid anchoring

Variable angle locking screw hole

Various screw insertion angles to secure rigid holding of ulnar head



Variable Angle Locking 30°

JEIL MEDICAL CORPORATION

702-703-704-705-706-804-805-807-812-ho, 55, Digital-ro 34-gil, Guro-gu, Seoul, 08378, Rep. of Korea

T 82.2.850.3500 F 82.2.850.3535 E jeil@jeilmed.co.kr

