

Mandible Reduction Dual Forceps


For Mandibular Bone Fracture Repair & Reconstruction

Craniomaxillofacial



Le Forte Mandible Reduction Dual Forceps

For Mandibular Bone Fracture Repair & Reconstruction




Stable & Safe Bone Holding for better surgical outcomes and shorter OP times

- Stable & strong holding force compared with conventional reduction forceps
- Wide joint motion range for various bone holding positions (up to 180° & 360° on each axis)
- Spike Tip for stable forceps purchasing on bone surface
- Improved holding stability with Ø2.0mm, 2.4mm & 2.7mm Screws
- Material: Biocompatible Stainless Steel

Parts List



Description	Part name
MRDF	111-158
T-Type Wrench	111-158-T
Flat Type Wrench	111-158-F



Screws hold the forcecp in fixed position

Optimal tip angle secures enough space to position mandibular reconstruction plate

JEIL MEDICAL CORPORATION

702-703-704-705-706-804-805-807-812-ho, 55, Digital-ro 34-gil, Guro-gu, Seoul, 08378, Rep. of Korea

DOMESTIC SALES_ T 82.2.850.3500 F 82.2.850.3535 E jeil@jeilmed.co.kr

INTERNATIONAL SALES_ T 82.2.850.3935 F 82.2.850.3537 E global@jeilmed.co.kr



LBEN-M9003(00)